



BLADE COLLECTION

CHART REFERENCES

LEGEND



Aggressive meniscal resection, fast removal of soft tissue for shoulder. Indicated for combination of soft tissue and bone resection.

Ø mm	CODE
3,5	FR1035LG
4,2	FR1042LG
4,8	FR1048LG
5,5	FR1055LG

NEPTUN



All purpose aggressive blade. Ideal for extensive tissue removal.

Ø mm	CODE
2,9	FR1029NP
3,5	FR1035NP
4,2	FR1042NP
4,8	FR1048NP
5,5	FR1055NP

STANDARD



Full radius: ideal for soft tissue resection and debridment.

Ø mm	CODE
3,5	FR1035ST
4,2	FR1042ST
4,8	FR1048ST
5,5	FR1055ST


EXPRESS



Oval bur. Fast bone resection. Indicated for subacromial decompression.

Ø mm	CODE
4,5	FR1045EX
6,0	FR1060EX

EPIC



Indicated for resistant tissues in knee.

Ø mm	CODE
3,5	FR1035EP
4,5	FR1045EP

ICON



Spherical bur. Ideal for notchplasty and osteochondral debridment.

Ø mm	CODE
4,5	FR1045IC
5,5	FR1055IC



ACCESSORIES COLLECTION

CHART REFERENCES

BATTERY STATION



CODE
FR20500

INTERNAL UNIT



CODE
FR20600

REMOVER



CODE
FR20700

INSERTER



CODE
FR20800

Contacts and worldwide distribution:

Frii Italy S.r.l. - Via Savona 97 - 20144 Milano

P +39 (0)2 422 969 05

E customerservice@friimedical.com - **W** www.friimedical.com