



THE FRII WAY *TO SURGERY*



THE 1st INTEGRATED **CORDLESS SINGLE-USE** SHAVER SYSTEM THAT MARKS A NEW GENERATION OF POWER INSTRUMENTS FOR ARTHROSCOPIC SURGERY.

In the 80's was created the first arthroscopic power instrument, since then all the arthroscopic companies tried to improve power instruments but never changed the concept.

Three decades later we believe it's time for a new perspective.

All the recent OR trends suggest

to reduce:

- size
- time
- maintenance

to improve:

- performance
- comfort
- simplicity

Our mission is to develop devices able to overcome the present standard, simplifying the use, reducing the maintenance, shortening the OR time, limiting the costs.



SUPERIOR FUNCTION BY DESIGN: THE ERGONOMIC AND CORDLESS HANDPIECE MINIMIZES USER FATIGUE AND OFFERS A MAXIMUM **FREEDOM** TO MANAGE ALL FUNCTIONS.

FEATURES & BENEFITS

CORDLESS AND WITHOUT CONSOLE

- all the functions are managed by hand control
- complete freedom to operate

STERILE SINGLE-USE DEVICE

- no preparation and no sterilization is required
- ready to use for saving time and money

REUSABLE INTERNAL UNIT

- high tech internal control unit
- complete rechargeable battery

ARTHROSCOPIC SURGERY FUNCTIONS

STANDARD FUNCTIONS

- forward and oscillating mode
- soft tissue and bony tissue resection
- adjustable fragments suction

ENHANCED FUNCTIONS

- seven pre-settings: to resume the ideal speeds for specific procedures
- light and ergonomic handpiece: superior comfort and more precision
- rotating blade: adjustable cutting-edge orientation on target



BLADE COLLECTION FRIIPOWER® SYSTEM OFFERS A WIDE ARRAY OF BLADES AND BURRS CAPABLE OF SUPPORTING EVEN THE MOST CHALLENGING ARTHROSCOPIC SURGICAL PROCEDURES.



ACCESSORIES COLLECTION FRIIPOWER® SYSTEM INCLUDES A SERIES OF ACCESSORIES THAT COMPLETE THE FUNCTIONS OF THE WHOLE EQUIPMENT.



All FRII® products are manufactured in certified facilities and are in pending to get compliance with the Food and Drug Administration (FDA) regulations. In addition, our products are CE marked according to the EEC Directives.


FRIIPOWER® IS THE BEGINNING OF A NEW APPROACH IN CAPITAL EQUIPMENT, MEETING THE CUTTING EDGE NEEDS IN TISSUE RESECTION.

FRII® is a new born company focused on arthroscopic surgery. The evolution of arthroscopic surgery has been one of the major developments in orthopedic surgery in the recent past: it has dramatically changed the orthopedic surgeon's approach to joint pathologies.

Following this path of innovation, FRII® has been created through extensive research and development programs conducted together with our partners and users as well as through a continuing refinement of existing technologies. Working together with surgeons, nurses and product specialists, FRII® is committed to develop high-quality products with increased benefits for patients and users as well as for their hospitals.

FRII® develops and manufactures medical devices with the mission to improve the quality of patient care and health care cost efficiency overcoming the present standard, paying a particular attention to comfort and ergonomic for the surgeon.





Authorized dealer

Contacts and worldwide distribution:

Frii Italy S.r.l. - Via Savona 97 - 20144 Milano

P +39 (0)2 422 969 05

E customerservice@friimedical.com - **W** www.friimedical.com



BLADE COLLECTION

CHART REFERENCES


LEGEND



| Ø mm | CODE |
|------|----------|
| 3,5 | FR1035LG |
| 4,2 | FR1042LG |
| 4,8 | FR1048LG |
| 5,5 | FR1055LG |

Aggressive meniscal resection, fast removal of soft tissue for shoulder. Indicated for combination of soft tissue and bone resection.


NEPTUN



| Ø mm | CODE |
|------|----------|
| 2,9 | FR1029NP |
| 3,5 | FR1035NP |
| 4,2 | FR1042NP |
| 4,8 | FR1048NP |
| 5,5 | FR1055NP |

All purpose aggressive blade. Ideal for extensive tissue removal.


STANDARD



| Ø mm | CODE |
|------|----------|
| 3,5 | FR1035ST |
| 4,2 | FR1042ST |
| 4,8 | FR1048ST |
| 5,5 | FR1055ST |

Full radius: ideal for soft tissue resection and debridement.


EXPRESS



| Ø mm | CODE |
|------|----------|
| 4,5 | FR1045EX |
| 6,0 | FR1060EX |

Oval bur. Fast bone resection. Indicated for subacromial decompression.


EPIC



| Ø mm | CODE |
|------|----------|
| 3,5 | FR1035EP |
| 4,5 | FR1045EP |

Indicated for resistant tissues in knee.

ICON



| Ø mm | CODE |
|------|----------|
| 4,5 | FR1045IC |
| 5,5 | FR1055IC |

Spherical bur. Ideal for notchplasty and osteochondrial debridement.



ACCESSORIES COLLECTION

CHART REFERENCES

BATTERY STATION



CODE

FR20500

INTERNAL UNIT



CODE

FR20600

REMOVER



CODE

FR20700

INSERTER



CODE

FR20800

Contacts and worldwide distribution:

Frii Italy S.r.l. - Via Savona 97 - 20144 Milano

P +39 (0)2 422 969 05

E customerservice@friimedical.com - **W** www.friimedical.com